

**Symphony
LabVIEW VI
Version 0.9.4
User Manual**

Preliminary

HORIBA JOBIN YVON

Table of Contents

Overview:.....	3
Minimum Requirements:	3
User_Lib.LLB	4
Spectral Mode.vi	4
Image Mode.vi	5
Version number.vi.....	6
Glob_Lib.LLB	7
Read Detector uniq-id.vi.....	7
Initialize CCD.vi	8
Init ADC & Gain.vi.....	10
Setup CCD.vi	11
Spectral Reformat CCD area.vi	14
Image Reformat CCD area.vi	15
Acquire image.vi.....	16
Acquire spectral.vi	18
Acquire.vi.....	20
Read Temperature.vi.....	21

Overview:

This set of LabVIEW VIs is designed to run a Symphony CCD system and provide the sub-VIs necessary to incorporate the Symphony CCD into a full LabVIEW application.

The LabVIEW library, USER_LIB.LLB contains two high level VIs to control a CCD and collect data. These are meant as examples for a user to follow and modify for using a Symphony CCD in their application.

GLOB_LIB.LLB contains the lower level sub-VIs that performs individual functions with the CCD. The individual VIs are used for the examples in USER_LIB.LLB and documented in the manual.

Minimum Requirements:

Windows 2000 or XP

LabVIEW 6.0.1 or Higher

Installed Symphony CCD system

Installed Symphony COM component – This can be installed either thorough Jobin Yvon's SynerJY software or the SynerJY SDK.

User_Lib.LLB

These are top level VIs that can be used as examples and a foundation for more advanced applications.

Spectral Mode.vi

This is an example program to show how to use the sub-VIs in an application. This can be used to collect a spectrum or a series of spectra from different CCD Areas.

Front Panel





Image Mode.vi

This is an example program to show how to use the sub-VIs in an application. This can be used to collect a spectrum or a series of spectra from different CCD Areas.

Front Panel



Version number.vi

This is a Top Level VI that holds Version information about this revision of the Symphony LabVIEW VIs.

Connector Pane



Front Panel




Glob_Lib.LLB

These are the low level sub-VIs that perform different functions with the Symphony CCD

Read Detector uniq-id.vi

This sub-VI will read all of the CCD information from the computer's registry. Registry settings are put into the registry on Hardware Configuration during installation of SynerJy software or the LabVIEW VIs.

Connector Pane



Front Panel



Controls and Indicators

[I32] **nbr of component** : [Output] Total number of Detector UniqueIDs found in the registry

[abc] **UniqueID** [Output] Array of UniqueIDs


[abc] **DeviceName** : [Output] Array of Device Names corresponding to the same index in Unique ID Array

[abc] **Name** : [Output] Array of Specific Names corresponding to the same index in Unique ID Array


Initialize CCD.vi

Sub VI to Initialize the CCD

Connector Pane



Front Panel



Controls and Indicators

- Unique id:** [Input] - This is the Unique ID of the CCD. This information can be found by using the Read Detector uniq-id.vi VI
- activeYPixels:** [Output] Total number of Y pixels of the CCD
- ActiveXPixels:** [Output] Total number of X pixels of the CCD
- Name:** [Output] Device Name of CCD
- Firmware Version:** [Output] Firmware version of the Symphony CCD Controller



Output error (Sortie d'erreur) The **error out** cluster can pass error information to other VIs in the software. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



status The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



code The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



source The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.




JYCCDLib.IJYCCDReqd


Link to the CCD COM component

Init ADC & Gain.vi

Connector Pane



Front Panel



Controls and Indicators



JYCCDLib.IJYCCDReqd



Input error (Entree d'erreur) The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



status The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



code The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



source The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



ADC description – [Output] Array of ADC Descriptions from the CCD Controller



ADC Number – [Output] Array of ADC Values from the CCD Controller




Gain Number – [Output] Array of Gain Descriptions from the CCD Controller

Setup CCD.vi

VI to setup a CCD for an acquisition

Connector Pane




Front Panel

The Front Panel of Setup CCD.vi includes the following controls and indicators:

- JYCCDLib.IJYCCDReqd**: A button to initiate the acquisition.
- DataSize**: A numeric control set to 0.
- Integration time**: A numeric control set to 10.
- Select ADC**: A dropdown menu set to JY_ADC_16_BIT.
- Gain**: A dropdown menu set to gain 1.
- Acquisition type**: A dropdown menu set to JYMCD_ACQ_FORMAT_SCAN.
- numAreas**: A numeric control set to 2.
- areas:** A table with 8 rows and 4 columns. The first column contains values: 0, 0, 1, 1024, 256, 1, 1, 0. The other three columns contain 0.
- entrée d'erreur (pas d'erreur)**: A status indicator showing a green checkmark and the code 00.
- sortie d'erreur**: A status indicator showing a green checkmark and the code 00.

Controls and Indicators

Gain : [Input] Set the Gain Level for the CCD Acquisition

 **Integration time** : [Input] Set the Integration Time for the Acquisition


 **Select ADC** : [Input] Select the ADC for the acquisition

 **numAreas** : [Input]


 **Acquisition type** [Input]

 **JYCCDLib.IJYCCDReqd**


This is a link to the CCD COM object functions

 **Input error (Entree d'erreur)** The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **status** The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.


The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.


 **code** The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.


 **source** The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.


 **areas** :

 **Output error (Sortie d'erreur)** The **error out** cluster can pass error information to other VIs in the software. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **status** The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **code** The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.




source The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

Spectral Reformat CCD area.vi

Connector Pane



Front Panel

areas :			
1	0	0	0
1	0	0	0
1	0	0	0
1024	0	0	0
256	0	0	0
1	0	0	0

Controls and Indicators:

numAreas [Input]

Total Number of Areas defined

X Pixels [Input]

Total Number of pixels in the X direction of the CCD

Y Pixels [Input]

Total Number of pixels in the Y direction of the CCD

Areas [Output]

Area Number: specific number of area for this set of parameters

X Start: First pixel in X dimension

Y Start: First pixel in y dimension

X Size : Total number of pixels in X direction

Y Size : Total number of pixels in Y direction


X binning : Number of pixels to sum together to get 1 data value in X direction

Y binning : Number of pixels to sum together to get 1 data value in Y direction.
For spectral information, this value needs to be equal to the Y size.

Image Reformat CCD area.vi

Sub-VI to define an area for acquiring an image from the CCD

Connector Pane



Front Panel

Parameter	Value
numAreas	1
X Pixels	1024
Y Pixels	256
Areas	0
Area number	0
X Start	1
Y Start	1
X Size	1024
Y Size	256
X Bining	10
Y Bining	10

Close

Controls and Indicators

numAreas [Input]

Total Number of Areas defined

X Pixels [Input]

Total Number of pixels in the X direction of the CCD

Y Pixels [Input]

Total Number of pixels in the Y direction of the CCD

Areas

Area Number: specific number of area for this set of parameters

X Start: First pixel in X dimension

Y Start: First pixel in y dimension

X Size : Total number of pixels in X direction

Y Size : Total number of pixels in Y direction


X binning : Number of pixels to sum together to get 1 data value in X direction

Y binning : Number of pixels to sum together to get 1 data value in Y direction


Acquire image.vi

This Sub-VI starts the acquisition of data from the Symphony CCD and returns a 3D array with X Data Value (Pixel), Y Data Value (Pixel), and Intensity

Connector Pane




Front Panel




Controls and Indicators

JYCCDLib.IJYCCDReqd


This is a link to the CCD COM object functions

 **Input error (Entree d'erreur)** The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **status** The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **code** The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **source** The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **ShutterOpen**

[I32]


True to open the shutter during acquisition
False to keep the shutter closed during an acquisition
type

data information from the image being acquired.


Acquire spectral.vi

This Sub-VI starts the acquisition of data from the Symphony CCD and returns a 2D array with X Data Value (Pixel) and Intensity

Connector Pane




Front Panel




Controls and Indicators

JYCCDLib.IJYCCDReqd


This is a link to the CCD COM object functions

 **Input error (Entree d'erreur)** The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **status** The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **code** The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **source** The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



ShutterOpen

True to open the shutter during acquisition

False to keep the shutter closed during an acquisition




type

data information from the image being acquired.


Acquire.vi

This Sub-VI starts the acquisition of data from the Symphony CCD and returns a 3D array with X Data Value (Pixel), Y Data Value (Pixel), and Intensity. This VI has the same function and Acquire Image.

Connector Pane




Front Panel




Controls and Indicators

JYCCDLib.IJYCCDReqd


This is a link to the CCD COM object functions

 **Input error (Entree d'erreur)** The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **status** The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **code** The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **source** The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.

 **ShutterOpen** : Input – Set to True to open the shutter during acquisition. Set to False to

keep shutter closed during acquisition.

I32


Type

Array of data reported back from detector. Data is in the format of X, Y, Intensity


Read Temperature.vi

Sub-VI to read the CCD's current temperature

Connector Pane



Front Panel



Controls and Indicators



JYCCDLib.IJYCCDReqd

This is a link to the CCD COM object functions



Input error (Entree d'erreur) The **error in** cluster can accept error information wired from VIs previously called. Use this information to decide if any functionality should be bypassed in the event of errors from other VIs.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



status The **status** boolean is either TRUE (X) for an error, or FALSE (checkmark) for no error or a warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



code The **code** input identifies the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



source The **source** string describes the origin of the error or warning.

The pop-up option **Explain Error** (or Explain Warning) gives more information about the error displayed.



units [Input]: Select Celsius, Fahrenheit, or Kelvin temperature units



Current Temperature [Output]: Shows the current CCD Temperature